

dima Print Teeth & Temp

Mechanical Properties of 3D-Printed Provisional Material Compared to Milled Material In vitro study – Kulzer R&D

3D-Printing technology and computer numerical control (CNC) milling machines are the most common computer-aided manufacturing methods in the dental field. While milling devices produce a workpiece using the subtractive technique, 3D-Printing is based on additive processes.¹

In addition to numerous other indications, fabricating long-term temporaries using computer-aided manufacturing methods is increasingly important. According to the corresponding requirements, these materials must exhibit stable and durable properties.²

This investigation compares the mechanical properties of a new photopolymer dima Print Teeth & Temp with milled provisional crown-and-bridge materials.

Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP

¹Park SM *et al.*: Flexural Strength of 3D-Printing Resin Materials for Provisional Fixed Dental Prostheses. Materials (Basel). 2020 Sep 8;13(18):3970. DOI: 10.3390/ma13183970.

²Gernet W *et al.*: Provisorische Versorgung. In: Zahnärztliche Prothetik. 2011. DOI: 10.1055/b-0034-18214

Mechanical Properties of 3D-Printed Provisional Material Compared to Milled Material

In vitro study – Kulzer R&D

Objectives

This investigation compares the mechanical properties of a new photopolymer with a milled provisional crown-and-bridge material by measuring flexural strength, elastic modulus, water sorption and solubility.


Methods

The tested provisional materials were dima Print Teeth & Temp (Kulzer GmbH), Telio CAD Disc and Telio CAD Block (both Ivoclar Vivadent GmbH).


Bar shaped specimens ($n=5$ /material) were produced for flexural strength testing and stored in water (24 h, 37 °C) before testing. 3-point bending was performed according to ISO 10477 with a universal testing device (Zwick/Roell Z010) at a cross-head speed of 1 mm/min to determine flexural strength. From the same measurement elastic modulus was calculated as secant (1N to 7N).

For water sorption and solubility testing, round specimens ($n=5$ /material) were produced from dima Print Teeth & Temp and Telio CAD Disc. Specimens were conditioned in a desiccator (22 h, 37 °C), placed in a second desiccator (2 h, 23 °C), weighed and procedure repeated until a constant mass was achieved. Specimens were weighed again after water storage (7 days, 37 °C), and after reconditioning. Testing and calculation according to ISO 10477. Statistics were done using one-way ANOVA (LSD $\alpha = 0.05$).

Results




Mean flexural strength (standard deviation) of Telio CAD Disc was 122.7(2.3) MPa, of dima Print Teeth & Temp was 113.8(2.3) MPa, and of Telio CAD Block was 103.5(5.8) MPa. The results were statistically significantly different between the materials (LSD $\alpha = 0.05$).



Mean elastic modulus (standard deviation) of dima Print Teeth & Temp was 3268(49) MPa, for Telio CAD Disc 3031(52) MPa and for Telio CAD Block 2538(87) MPa. The results were statistically significantly different between the materials (LSD $\alpha = 0.05$). For elastic modulus no ISO requirement is stated.

Mechanical Properties of 3D-Printed Provisional Material Compared to Milled Material

In vitro study – Kulzer R&D



Water sorption and solubility were not measured for Telio CAD Block. The results of water sorption were statistically significantly different between the materials (LSD $\alpha=0.05$). All results of water sorption and solubility fulfil the ISO 10477 requirements.

Conclusion

Within the limitations of this study, it can be stated that the new photopolymer dima Print Teeth & Temp exhibits good mechanical properties comparable to those of a milled reference product.

Comment

Restorations are not only exposed to several forces in the mouth. They are also subjected to the oral environment's conditions, influencing materials' properties and longevity.³ dima Print Teeth & Temp presents excellent resistance to fracture and elastic deformation due to its high flexural strength and elastic modulus value.⁴ Water sorption may affect composite materials by reducing the wear resistance; for this reason, it is necessary that these materials uptake the least possible amount of water. Low solubility is also essential, meaning less mechanical degradation and release of residual monomers.⁵

Source

Derendorf V, Kastrati A, Loh W, Herr A: Poster 278 at the PER-IADR Oral Health Research Congress in Marseille, September 15–17, 2022.

The study was abbreviated, summarised and commented and all diagrams and titles have been established by Kulzer.

³Sideridou ID: Aging studies of light cured dimethacrylate-based dental resins and a resin composite in water or ethanol/water. Dental Materials. 2007;23:1142-9

⁴Arnetz G: Neue Materialien, neue Perspektiven. Zahntech Mag 2015; 19(1): 18-25

⁵Oysaed H: Composites for use in posterior teeth: mechanical properties tested under dry and wet conditions. J Biomed Mater Res. 1986 Feb; 20(2): 261-71. doi: 10.1002/jbm.820200214