Diamond Burs

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GB Diamond Burs

Rotary instruments for dental application

Safety recommendation

Handpieces must be in impeccable technical and hygienic condition. Ensure the cooling jet nozzles are unlocked and free from debris at all times. The FG burs (instruments) must be inserted into the handpiece carefully, without the use of force. They must be inserted axially and as deep as possible. Discard blunt, deformed or damaged instruments immediately. Improper use leads to increased risk as well as greater wear, with potentially substandard results. Please observe maximum permissible speeds. To avoid unnecessary damage to the substrate being worked on, ensure the bur is running at the recommended operating speed **before** cutting begins.

Note

Using burs with longer shafts or thin necks (e.g. microdissection instruments) may have a tendency to resonate (vibrate) or deflect when the maximum permitted rotational speed is exceeded. This can lead to breakage or severe damage to the substrate. Non-observance will result in an increased safety risk! Therefore, it is strongly recommended to observe the maximum speed on the package labels.

Applying force

The forces applied by the operator, from the bur to the substrate, should be low and not exceed 2 N maximum. Excessive force can lead to breakage of the diamond coating or deflection of the instrument, as well as an increase in heat generation between the bur and the substrate. Tilting and leverage of the instrument must be avoided as it leads to increased risk of breakage. **Excessive force can lead to undesirably rough surfaces and/or damage to the pulp as a result of overheating. In extreme cases, instrument breakage may also occur.**

Cooling

To avoid unnecessary damage to the natural tooth substance as well as restorations, continuous and sufficient **water spray cooling** must be ensured throughout use. In the speed range >1,500 min-1, a minimum of 50 ml/min is necessary to ensure reliable cooling of the bur and substrate. The cooling spray facility of the handpieces are tailored to ISO standard sizes. Therefore, additional cooling is necessary for FG instruments with an overall length >19 mm and with a working part diameter >2 mm.

If there is insufficient water cooling, this may lead to irreversible damage of the teeth and surrounding tissue.

Maximum speeds: FG instruments (ISO 314)

Diameter of the working part ISO classification	Maximum permissible speed (min-1)
005 - 010	450.000
012 - 016	450.000
018 - 021	300.000
023 - 031	160.000
033 - 055	120.000

Maximum speeds: Handpiece instruments (ISO 104/204)

Diameter of the working part ISO classification	Maximum permissible speed (min-1)
003 - 023	50.000
025 - 040	40.000
040 - 060	30.000
060 - 070	25.000
070 - 250	25.000
300 - 380	20.000
450	15.000

Optimum speed for all instruments usually 40 to 50 % of the maximum specification.

Please observe the information on the package labels of the instruments. Exceeding the specified maximum speed can lead to a reduction in performance of the instruments and/or could endanger the patient or operator.

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Reprocessing

Manufacturer information regarding reprocessing in accordance with the RKI guideline “Infection Prevention in Dentistry, Hygiene Requirements“ Medical Devices Semi-Critical B and Critical B according to DIN EN ISO 17664.

Products

This manufacturer information applies to all diamond cutting instruments supplied by Kulzer which are used for surgical, periodontal or endodontic procedures. These are rotary diamond instruments. The product range only comprises of instruments that have been supplied un-sterile. These must be processed prior to their first and all subsequent uses.

Limitation of reprocessing

The end of a product’s lifetime is usually determined by wear and damage as a result of use. Frequent reprocessing has no effect on the performance of these instruments. Disposable products (marked with  on the label) must not be re-used or reprocessed.

Storage and transport

Immediately after use, or at the latest 1 hour after use on the patient, insert instruments into a suitable bur disinfector filled with cleaning agent/disinfectant (e.g. BIB forte eco). Bur disinfectant is manufactured in a concentration according to the manufacturer’s instructions, with e.g. BIB forte eco: prepare the concentrate with water, first add water, cover bur disinfector. Observe reaction time (e.g. BIB forte 0.5 % 60 min). The instruments should be transported to the processing location, protected from contamination, in the bur disinfector.

Cleaning and disinfection

In accordance with the recommendation of the Robert Koch Institute (RKI), Semi-Critical B Products should preferably be processed mechanically. Critical B Products should generally be processed mechanically. In the case of products with long, narrow lumina or non-smooth surfaces, cleaning must take place mechanically. Silicone stoppers must be removed before the processing of instruments.

Mechanical cleaning

Equipment

Cleaning brushes, cleaning/disinfection device (RDG) (e.g. Miele with Vario TD programme), in accordance with EN ISO 15883, cleaning agent (e.g. 0.5 % Neodisher mediclean cleaner), instrument stands for rotary, oscillating instruments (stainless steel instrument stands).

Procedure

Take instruments out of the instrument stands/bur disinfector immediately before mechanical processing and brush within the cold water bath until all visible contamination has been removed. Load the instruments into the instrument stands. Start mechanical cleaning while observing the manufacturer’s instructions for the RDG and cleaner. The following process is validated: Vario TD programme: 2 min. pre-cleaning, 5 min. cleaning at 55°C with a cleaning agent, 3 min. neutralisation, 2 min. flushing, final rinse with demineralised water 5 min. at >90°C.

Thermal disinfection in the RDG

Carry out mechanical cleaning of the products fixed in the instrument stands/bur disinfector (e.g. RDG Miele with Vario TD programme) incl. thermal disinfection. Adhere to the manufacturer’s instructions. With a validated RDG, disinfection is verifiably guaranteed. Kulzer products are thermostable to 134°C.

Manual cleaning and disinfection – standardised procedure

Equipment

Cleaning brushes (e.g. plastic brushes, sterilisable), ultrasonic cleaner, cleaning agent and disinfectant for dental instruments with tested efficacy (e.g. BIB forte eco, Alpro Medical), instrument stands for rotary, oscillating instruments (e.g. stainless steel instrument stands). Adhere to the manufacturer’s instructions.

Procedure

Take instruments out of the instrument stands/bur disinfector immediately before manual processing and brush within the cold water bath until all visible contamination has been removed. Rinse instrument and instrument stand under running water. Place instruments into a suitable sieve container within the ultrasonic device filled with cleaning agent and disinfectant. Carry out cleaning and disinfection in accordance with the instructions of the manufacturer of the ultrasonic cleaner and the cleaning agent and disinfectant; e.g. BIB forte eco 3 % - 10 min. (at 55°C tested according to EN 14476). Rinse instrument thoroughly with suitable water (e.g. demineralised water) once the reaction time has elapsed. Preferably dry instruments with medical compressed air. In accordance with KRINKO, manual processing is completed by means of thermal disinfection in the steam steriliser. Adhere to the manufacturer’s instructions.

Visual inspection for integrity and cleanliness with an appropriate magnification device (8-10-fold magnification is recommended). If residual contamination is still detectable on the instrument after processing, repeat cleaning and disinfection until contamination is no longer visible. Instruments which exhibit defects must be disposed of immediately, e.g. missing diamond coating, blunt and flared cutting, distorted shapes, corroded surfaces or residual contamination which cannot be removed.

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Stand: 2017-11

Final processing steps

Medical products Critical B: Sterile packaging and sterilisation

Packaging

Prior to sterilisation, the instruments must be placed within the instrument stands and together packed in two transparent sterilisation bags (e.g. Steriking or VP Steriking) and heat sealed with the seam sealer (e.g. Hawo). The instruments must be protected. A suitable standardised procedure must be used for packaging.

Sterilisation

Successful sterilisation of the packaged instruments is successfully demonstrated in the prevacuum-steam sterilisation procedure (e.g. steam steriliser, MMM Selectomat HP) with the following minimum parameters: 3 prevacuum phases, 132°C sterilisation temperature, holding time 3 min. (full cycle), drying time 10 min. The equipment manufacturer’s instructions must be observed. **Note: The products are not suitable for sterilisation in the chemiclave and hot air disinfectior.**

Transport and storage

The processed instruments should be transported and stored protected from recontamination. In the case of sterile materials, they should also be protected from dust and moisture.

Basics

Observe the applicable legal provisions for the reprocessing of medical products in your country (e.g. www.rki.de). The manufacturer guarantees that the aforementioned processing procedures for the processing of the aforementioned instrument groups are appropriate for their re-use according to their intended purpose. The operator is responsible based on his risk assessment, for the actual reprocessing performed achieving the specified objectives for the intended use with the equipment, materials, process parameters and personnel employed.

Routine inspections of the validated mechanical or the standardised manual processing procedures are necessary for this. Likewise, any deviation from the validated procedures should be carefully evaluated and approved by the processor in terms of efficacy and potential adverse consequences.

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