

New twist  
shapes



## Venus® supra

### The simple way to high gloss

Long-lasting aesthetics are what Venus does best. With the polishing system Venus supra dentists add a wonderful high gloss to the finished composite restorations.

Venus supra now comes with new benefits:

- » **Two new twist shapes:** The Twist Disc is ideal for polishing results on various composite surfaces. The Twist Brush achieves perfect polishing results on e.g. occlusal surfaces with easy access to deep fissures.
- » **New separate refills:** Pre – and high-gloss polishers can now be ordered separately. This helps optimize stock-keeping and re-order cycles.

Giving a hand to oral health.



**KULZER**  
MITSUI CHEMICALS GROUP

# Venus® Supra – the long-lasting high-gloss polisher system:

- Two-step polishing system consists of pre- and high-gloss variants in different shapes
- Made of silicone highly filled with microfine diamond powder
- Tailor-made for Venus composites such as Venus Diamond ONE and Venus Pearl ONE
- Universal solution for all anterior and posterior restorations, giving them an aesthetic shine and smooth surfaces
- Autoclavable and reusable

## New twist shapes

- Ease polishing and adaptation to natural tooth morphology

## New refills


- Optimise stock keeping and reorder cycles as pre and high-gloss polishers can be ordered separately




Ulf Krueger-Janson

*"I've been using Venus Supra polishers for over 13 years. They are indestructible and enable the best high gloss that I have ever been able to achieve. Simply great!"*

### Venus Supra polishing leads to excellent polishing results




Versatility: Twist Disc – ideal polishing results on various composite surfaces



Art. Code	Article	Content
<b>Intro Kit</b>		
66087541	Venus Supra Intro Kit 2x5	1 pre and 1 high gloss polisher of each shape (Flame, Cup, Large Flame, Twist Brush, Twist Disc)
<b>New Shapes</b>		
66093899	Venus Supra Refill Pre Twist Disc	10 pre polishers of shape Twist Disc
66093900	Venus Supra Refill Gloss Twist Disc	10 high gloss polishers of shape Twist Disc
66093901	Venus Supra Refill Pre Twist Brush	10 pre polishers of shape Twist Brush
66093902	Venus Supra Refill Gloss Twist Brush	10 high gloss polishers of shape Twish Brush
<b>New Refills</b>		
66093893	Venus Supra Refill Pre Cup	10 pre polishers of shape Cup
66093894	Venus Supra Refill Gloss Cup	10 high gloss polishers of shape Cup
66093895	Venus Supra Refill Pre Flame	10 pre polishers of shape Flame
66093896	Venus Supra Refill Gloss Flame	10 high gloss polishers of shape Flame
66093897	Venus Supra Refill Pre Large Flame	10 pre polishers of shape Large Flame
66093898	Venus Supra Refill Gloss Large Flame	10 high gloss polishers of shape Large Flame



Ideal results on fissures: Twist Brush



The brush bristles can be bent for easy access to fissures

[kulzer.com/try-more-venus](http://kulzer.com/try-more-venus)