

artic 

GB Mould chart
DE Formenkarte
FR Cartes des Formes
ES Carta de Formas
IT Tabella forme
PT Carta Molde
NL Vormenkaart
SE Formkort
DK Formkort
NO Formkort
FI Muotto korti
LV Formas karte
LT Formos kortelė
PL Karta kształtu
HR Kartica oblika
RU Карточка формы
SK Tvarová karta
BA Kartica oblika
CZ Karta tvaru
SI Kartica oblike
BGR плесен диаграма
ARA مخطط القالب



Artic®

Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP

artic 6 Anteriors

Frontzähne/ Antérieures / Anteriores / Anteriori / Anteriores / Tanden / Framtänder / For tænder / Fortennene / Etuhampaat / Priekšējie zobi / Priekiniai dantys / Zęby przednie / Prednji zubi / Передние зубы / Predné zuby / Predni zubi / Přední zuby / Sprednji zobje / предни / الأمامية

O 10 	40,9 9,0 7,5	S 35 	50,6 10,8 9,2
O 20 	46,1 9,7 8,4	S 40 	53,2 12,6 9,9
O 30 	47,0 11,0 8,6	T 10 	44,1 9,0 8,2
R 20 	43,0 10,0 8,0	T 15 	42,6 9,6 7,9
R 30 	47,8 11,8 8,9	T 20 	45,6 9,8 8,8
R 40 	51,0 13,4 9,7	T 30 	45,5 10,6 8,2
S 05 	40,3 8,9 7,4	T 35 	47,4 10,0 8,8
S 10 	44,2 9,1 7,9	T 40 	49,1 12,1 8,7
S 20 	44,0 9,9 7,9	TO 05 	42,1 9,0 7,8
S 30 	49,8 10,5 8,9	TO 10 	43,3 9,3 8,3

artic 6 Anteriors

Frontzähne/ Antérieures / Anteriores / Anteriori / Anteriores / Tanden / Framtänder / For tænder / Fortennene / Etuhampaat / Priekšējie zobi / Priekiniai dantys / Zęby przednie / Prednji zubi / Передние зубы / Predné zuby / Predni zubi / Přední zuby / Sprednji zobje / предни / الأمامية















artic 8 Posteriors 10°/20°





























Seitenzähne / Postérieures / Posteriores / Posteriori / Posteriores / Kiezen / Bakre tänder / Bageste tænder / Bakre tenner / Takimmaisiet hampaat / Aizmugurējie zobi / Užpakaliniai dantys / Zęby tylne / Stražnji zubi / Задние зубы / Zadné zuby / Stražnji zubi / Zadní zuby / Zadnji zobje / постериори / الخلفية


TO 20		43,8 9,6 7,9	L 25		37,0 9,2 5,7
TO 30		45,9 10,4 8,6	L 30		35,9 9,0 5,5
TO 40		47,8 11,0 8,9	L 40		36,0 9,0 5,5
L 05		31,5 9,0 5,0	L 45		38,9 9,6 5,9
L 10		33,5 7,8 5,1	L 50		39,2 11,1 5,9
L 15		33,9 9,5 5,2	L 55		40,1 10,6 6,1
L 20		35,5 8,9 5,3			
30 U	0°			28,9 8,4	
30 L				30,7 9,1	
32 U	0°			33,1 8,8	
32 L				34,8 9,2	

artic 8 Posteriors 10° / 20°

Seitenzähne / Postérieures / Posteriores / Posteriori / Posteriores / Kiezen / Bakre tänder / Bageste tænder / Bakre tenner / Takimmaisest hampaat / Aizmugurējie zobi / Užpakaliniai dantys / Zęby tylne / Stražnji zubi / Задние зубы / Zadné zuby / Stražnji zubi / Zadní zuby / Zadnji zobje / постериори / الخلفية

30 U	10°			29,4
30 L				8,8
32 U	10°			32,1
32 L				9,5
30 U	20°			29,2
30 L				7,6
32 U	20°			31,0
32 L				6,9
34 U	20°			32,5
34 L				9,1
36 U	20°			35,3
36 L				8,2
36 U	20°			34,0
36 L				10,0
				37,6
				8,9

	Upper anteriors	Lower anteriors	Posteriors	0°	10°	20°
	 O 10	L 10		30	30	30
	 O 20	L 30		30, 32	30, 32	32
	 O 30	L 40		32	32	32
	 R 20	L 20		32	32	32, 34
	 R 30	L 25		32	32	36
	 R 40	L 55		32	32	36
	 S 05	L 05		30	30	30
	 S 10	L 10		30	30	30
	 S 20	L 20		30, 32	30, 32	32
	 S 30	L 40		32	32	32, 34
	 S 35	L 55		32	32	36
	 S 40	L 55		32	32	36
		 T 10	L 10		30	30
 T 15		L 15		30, 32	30, 32	32
 T 20		L 30		30, 32	30, 32	32
 T 30		L 20		32	32	32
 T 35		L 40		32	32	32
 T 40		L 50		32	32	34
	 TO 05	L 05		30	30	30
	 TO 10	L 15		30	30	30
	 TO 20	L 20		30, 32	30, 32	32
	 TO 30	L 40		32	32	32
	 TO 40	L 45		32	32	34



This mould chart complies with EN ISO 22112:2017. / Diese Formenkarte erfüllt EN ISO 22112:2017. / Cette carte forme remplit EN ISO 22112:2017. / Esta carta de formas es conforme con la EN ISO 22112:2017. / Questa tabella forme soddisfa EN ISO 22112:2017. / Esta carta molde está em conformidade com EN ISO 22112:2017. / Deze vormenkaart voldoet aan EN ISO 22112:2017. / Detta formkort uppfyller EN ISO 22112:2017. / Dette formkort opfylder EN ISO 22112:2017. / Dette formkortet oppfyller EN ISO 22112:2017. / Tämä muotokortti täyttää EN ISO 22112:2017. / Ši formos korta atbilst EN ISO 22112:2017. / Ši formos korta išpildo EN ISO 22112:2017. / Ta karta ksztaftu spełnia EN ISO 22112:2017. / Ova kartica oblika ispunjava EN ISO 22112:2017. / Эта карта формы выполняет EN ISO 22112:2017. / Táto tvarová karta spĺňa EN ISO 22112:2017. / Ova kartica oblika ispunjava EN ISO 22112:2017. / Tato karta tvaru splňuje EN ISO 22112:2017. / Ta oblika kartice izpolnjuje EN ISO 22112:2017. / Тази диаграма на формата е в съответствие с EN ISO 22112:2017. / يتوافق مخطط القوالب هذا مع EN ISO 22112:2017.

**Manufacturer**

Kulzer GmbH
Leipziger Straße 2
63450 Hanau
Germany

Contact in Germany

Kulzer GmbH
Leipziger Straße 2
63450 Hanau
Germany
info.lab@kulzer-dental.com