

DIMA PRINT TEETH & TEMP**2-IN-1: NATURAL LOOK &
FAMILIAR HANDLING****NEW**DISCOVER
MOREKULZER.COM/
TEETHTEMP

With **dima Print Teeth & Temp** you can print two indications at the same time with your cara printer: teeth & temporaries and you get the same benefits for both. dima Print Teeth & Temp delivers excellent properties, convincing benefits and a natural look for both indications.

**KULZER**
MITSUI CHEMICALS GROUP

DIMA PRINT TEETH & TEMP**2-IN-1: NATURAL LOOK & FAMILIAR HANDLING****Benefits:**

- ▶ dima Print Teeth & Temp material covers two applications at once: temporaries and teeth for removable dentures. Less material changes, less cleaning, less storage.
- ▶ Approx. 25 min. to print teeth and approx. 50 min. to print temps (or temps + teeth together¹).

SAVE TIME AND MONEY

- ▶ No extra mile! The handling of the finishing process is comparable with conventional production methods, like milling.

CONVENIENT**EFFICIENT PRINTING**

- ▶ Excellent mechanical properties. A study recreating the conditions in the oral cavity showed that dima Print Teeth & Temp is a reliable material. Learn more at: kulzer.com/teethtemp-science

COMPARABLE TO MILLED MATERIALS

¹ Printing time depends directly on the resolution chosen. Temps require a higher resolution than teeth.



Human tooth



Kulzer natural look










Watch Video –
Natural look
[kulzer.com/
teethtemp-natural-look](http://kulzer.com/teethtemp-natural-look)

Features:

- ▶ Four color blends developed to cover 7 common shades*
- ▶ Discover natural fluorescence
- ▶ Very low minimum wall thickness of 0.3 mm possible

dima Print Teeth & Temp shades

Product	7 shades
dima Print Teeth & Temp A1/B1	 A1
	 B1
dima Print Teeth & Temp A2/B2	 A2
	 B2
dima Print Teeth & Temp A3/A3.5	 A3
	 A3.5
dima Print Teeth & Temp BL	 BL

*Color may vary depending on the used equipment.

Overview wall thickness**Note:**

For more details, please contact your local Kulzer representative.

