

Aligner Attachments and Retainer Fixation

No switching. No compromise. Just seamless integration.



I'm placing aligner attachments regularly, but ...

"I don't want to switch materials just for attachments."

"I need a reliable way to place attachments and remove them after treatment without damaging enamel."

"I want to keep my workflow simple and efficient."

Venus® composites are now officially approved for aligner attachments and retainer fixation.

This allows you to keep things simple: use the material you trust for both restorative and orthodontic procedures, without the hassle of switching products or established routines. Whether you prefer Venus Bulk Flow ONE, Venus Pearl, Venus Diamond or Venus Diamond Flow – you can now rely on your go-to composite for even more indications.

Why it matters:

- **Streamlined inventory**
One material for restorative and orthodontic treatments = less stock and more simplicity.
- **Familiar handling**
Save time by using your trusted composite, no additional training needed.
- **ONE shading**
Enjoy the convenience of a single shade that blends seamlessly with any tooth colour. No need to change shades, even for darker canines.
- **Confident removal**
The Kulzer Translux Cure in DETECT mode helps detect any residues thanks to the intrinsic fluorescence of Venus composites – for clean, confident results.

Choose it. Use it. Love it.



KULZER
MITSUI CHEMICALS GROUP

Aligner attachments by Dr Amisha Patel, UK



Aligner buttons made of Venus Pearl ONE shade.



Aligner in situ.



"The ONE shade helps keep the process of straightening teeth with clear aligners discreet. It eliminates the need for shade matching and offers every patient 'barely there buttons'."

Dr Amisha Patel, UK

Get started right away with the new
Venus Bulk Flow ONE Aligner Attachment Kit.

Art. Nr. 66112380

Contains:

- 1x 2g syringe Venus Bulk Flow ONE
- 1x 2ml bottle iBOND Universal
- 5 cannulas



Discover more:



kulzer.com/venus



Workflow tip

For secure attachment placement, apply **iBOND® Universal** in self-etch mode – ensuring strong adhesion while protecting enamel.

After treatment, polish with **Venus Supra** to gently remove any residual composite and achieve a smooth, natural finish.

Smart placement. Gentle removal. Seamless results.

